



## Large Classic BPA Free Silicone Hot Water Bag For Neck Shoulder And Hand Feet Warmer

### Our Product Introduction

for more products please visit us on [siliconeproducts-factory.com](http://siliconeproducts-factory.com)

#### Basic Information

- Place of Origin: China
- Brand Name: JY
- Certification: FDA, LFGB
- Model Number: JY-255
- Minimum Order Quantity: 50pcs
- Price: negotiable
- Packaging Details: opp bag, customized pakaging
- Delivery Time: negotiable
- Payment Terms: T/T
- Supply Ability: negotiable



#### Product Specification

- Material Composition: Food Grade Silicone
- Tasteless: Yes
- Color: Pink, Blue, Red
- Design Style: Modern
- Brand: JY
- Keywords: Silicone
- OEM: Accepted
- Highlight: **BPA Free Silicone Hot Water Bag, Silicone Hot Water Bag For Neck**



#### More Images



## Product Description

# Hot Water Bottle with Hand Pocket Cover

### Product description

Our Hot Water Bottle is made from advanced thermoplastic material, ensuring it is long-lasting and sturdy. Treats Neck and Shoulder Pain Relief, Headaches, Cramps, Back Pain, Sore Muscles, Injuries and more! Its versatility makes it an essential item for all seasons.

### Product specification

Special Feature	Color Changing
Application	Hot
Size	S, M, L
Tasteless	Yes
BPA free	Yes
Color	Customized

### Product Features

- [Professional Material] -- Thermoplastic material which is odorless and recyclable. This kind of thermoplastic material can hold heat longer than a traditional hot water bottle. Could be used as an ice pack, or cold water bottle. Note, this is not suitable for Microwave oven heating.
- [Safe & Easy to Use] -- The hot water bag is explosion-proof, leak-proof which is good for daily use. Pregnant women and children can use it safely since it is odorless. Please fill only two-thirds 80 degree water into it.
- [Product Highlights] -- The hot water bottle has been done the leak detection and sealing detection within 24 hours before leaving the factory. Also, the product can retain heat for an extra long time. Ribbed outer surface helps maintain liquid temperature.
- [Best Assistant] -- This hot water bottle can ease joint disease, sore muscles, cold and cough, dysmenorrhea and other common aches. It can also be a cold compress for people with fever or to reduce swelling for minor injuries.

### Product More Details

- [Premium Material] -- The product is made of high quality PVC material, non-toxic and odorless. Compared to the traditional rubber hot water bags, this one is more durable and longer lasting insulation. Note: Please do not use microwave heating.
- [Safe & Easy to Use] -- Hot water bags have passed 24 hours of explosion-proof and leak-proof test before leaving factory. BPA free material is more friendly to pregnant women, the elderly and children. Please fill two-thirds of 80 degrees of water into it.
- [Unique Design] -- The special pattern design on the surface of the hot water bag, high-quality raw materials, and exquisite craftsmanship, making the hot water bag has a more durable insulation effect.
- [Multi Functional] -- Hot water bag can be used as hot or cold compress for relieving muscle aches and pains, relieving menstrual cramps, sports injuries, etc. • Long-lasting Heat Retention: Unlike traditional hot water bottles, our innovative thermoplastic material provides superior insulation, allowing the hot water to stay warm for a longer period of time. Say goodbye to constantly reheating your hot water bottle and enjoy consistent heat therapy for hours.
- Secure and Leak-proof Design: The Hot Water Bottle features a secure seal and a leak-proof cap, ensuring that the hot water stays contained and prevents any accidental spills or burns. You can confidently use it as a heat pack without worries.
- Multi-purpose Usage: It can be used as a cold pack or an ice bag for soothing relief from swelling, sports injuries, or general cooling therapy. Its versatility makes it an all-in-one solution for various pain management needs.
- Caution: Made Of Natural Rubber Which Smells; To Reduce The Smell of Rubber Use Vinegar/Water Mix and Submerge For Min of 30min











**Silicone JinYu Industrial Co., Ltd.**

 +86 15028563200

 [angela@siliconeproducts-factory.com](mailto:angela@siliconeproducts-factory.com)

 [siliconeproducts-factory.com](http://siliconeproducts-factory.com)

Room 306, No. 3 Shengyuan Street, Yayuan, Nancheng Street, Dongguan China